


Uninstall existing RadStar Data Capture Program


Note: Windows Vista Instructions where different are shown [xxxxx "XXXXX"]

1. Click "**Start**" and click "**Control Panel**".

Windows XP




Vista




2. Double click ***"Add/Remove Programs"***.

Windows XP




[Click "*Uninstall a program*".]

Vista




- Click **RadStarDataCapture** and click **Remove** button.

Windows XP




[Scroll down, click “RadStarDataCapture” and click “Uninstall” button].

Vista




4. A window will appear asking if you are sure you want remove RadStarDataCapture, click **"Yes"**.

Windows XP



[A window will appear asking if you are sure you want to uninstall RadStarDataCapture, click "Yes". Click "Allow" in the pop up User Account Control Window (not shown)].

Vista




5. Wait for Windows to complete Uninstalling and then close "Add or Remove Programs" and "Control Panel" windows.


SUPPLEMENT 1 – WINDOWS 7 INSTRUCTIONS

Uninstall existing RadStar Data Capture Program


1. Display the Windows **"Start"** menu by clicking on Windows "Start" icon  .
2. Click on the **"Control Panel"** item on left side of "Start" menu.




3. On **“Control Panel”** menu click on the **“Programs and Features”** menu item.



4. On the “Programs and Features” menu double click on the “**RadstarDataCapture for Windows Vista/7**”.



5. If the confirmation window shown below appears click on **“Yes”**.



6. After the uninstall process completes the **“RadStarDataCapture for Windows Vista/7”** software will be uninstalled and the program name will be deleted from the “Programs and Features” list.